




PARAMCARE LIFESCIENCES PVT. LTD.
where innovation starts...

RapiChex[®]

PCP RAPID TEST KIT



FEATURES:

- ✓ Easy to use with all accessories.
- ✓ Specimen type 2 drops (50µL) urine samples.
- ✓ Results within 5 mins & it works on reverse principle.
- ✓ Product Presentation: 10/25/50 Tests.
- ✓ Sensitivity: 100%, Specificity: 98.57%.
- ✓ Room Temperature Storage.
- ✓ Cut off: 10 ng/ml.



Contact@paramcarelifesciences.com
Directors@paramcarelifesciences.com



WWW.PARAMCARELIFESCIENCES.COM